



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

COPS5 Antibody [#C11442]

Catalog Number: C11442

Concentration: 1mg/ml

Swiss-Prot No.: Q92905

Other Names: EC 3.4.-.-, Signalosome subunit 5, SGN5, Jun activation domain-binding protein 1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human COPS5.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: COPS5 antibody detects endogenous levels of total COPS5 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Claret F.-X., Nature 383:453-457(1996).

Asano K., J. Biol. Chem. 272:27042-27052(1997).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

