



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## TAF1A Antibody [#C11088]

**Catalog Number:** C11088

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q15573

**Other Names:** TATA box-binding protein-associated factor 1A, TBP-associated factor 1A, TBP-associated factor RNA polymerase I 48 kDa, TAFI48, Transcription factor SL1

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human

TAF1A.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** TAF1A antibody detects endogenous levels of total TAF1A protein.

**Reactivity:** Human, Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:5000

**References:**

Comai L., Science 266:1966-1972(1994).

Ota T., Nat. Genet. 36:40-45(2004).

Gregory S.G., Nature 441:315-321(2006).

