



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TSN Antibody [#C11090]

Catalog Number: C11090

Concentration: 1mg/ml

Swiss-Prot No.: Q15631

Other Names: Translin, TSN

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human TSN.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TSN antibody detects endogenous levels of total TSN protein.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 IHC: 1:50~1:100
ELISA: 1:40000

References:

Aoki K., Nat. Genet. 10:167-174(1995).


Hillier L.W., Nature 434:724-731(2005).

Aoki K., Genomics 43:237-241(1997).



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human brain tissue using TSN antibody.



Peptide - +

Western blot analysis of extracts from Jurkat cells, using TSN antibody.