



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TRIP4 Antibody [#C11094]

Catalog Number: C11094

Concentration: 1mg/ml

Swiss-Prot No.: Q15650

Other Names: Activating signal cointegrator 1, ASC-1, Thyroid receptor-interacting protein 4, TRIP-4, TRIP4

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human TRIP4.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TRIP4 antibody detects endogenous levels of total TRIP4 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Kim H.J., Mol. Cell. Biol. 19:6323-6332(1999).

Lee J.W., Mol. Endocrinol. 9:243-254(1995).

Jung D.-J., Mol. Cell. Biol. 22:5203-5211(2002).

