



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## FOXJ1 Antibody [#C11446]

**Catalog Number:** C11446

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q92949

**Other Names:** Forkhead-related protein FKHL13,  
Hepatocyte nuclear factor 3 forkhead homolog 4, HFH-4

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human FOXJ1.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** FOXJ1 antibody detects endogenous levels of total FOXJ1 protein.

**Reactivity:** Human, Mouse, Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Pelletier G.J., Am. J. Physiol. 274:L351-L359(1998).

Murphy D.B., Genomics 40:462-469(1997).

Li C.-S., J. Hum. Genet. 51:292-297(2006).

