



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## DLX4 Antibody [#C11450]

**Catalog Number:** C11450

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q92988

**Other Names:** Homeobox protein DLX-4, DLX-7, DLX-8, Beta protein 1

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from C-terminal of human DLX4.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** DLX4 antibody detects endogenous levels of total DLX4 protein.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:1000

**References:**

Nakamura S., Genomics 38:314-324(1996).

Price J.A., Hum. Mol. Genet. 7:563-569(1998).

Haga S.B., Leukemia 14:1867-1875(2000).

