



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

NAB2 Antibody [#C11102]

Catalog Number: C11102

Concentration: 1mg/ml

Swiss-Prot No.: Q15742

Other Names: EGR-1-binding protein 2; MADER;
Melanoma- associated delayed early response protein;
NGFI-A-binding protein 2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human NAB2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: NAB2 antibody detects endogenous levels of total NAB2 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 IHC: 1:50~1:100


ELISA: 1:10000

References:

Svaren J., Mol. Cell. Biol. 16:3545-3553(1996).

Kirsch K.H., Oncogene 12:963-971(1996).


Svaren J., Genomics 41:33-39(1997).



Peptide -

+

Immunohistochemistry analysis of paraffin-embedded human brain tissue using NAB2 antibody.



Western blot analysis of extracts from HT-29 cells and RAW264.7 cells, using NAB2 antibody.