



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

SPR1 Antibody [#C11103]

Catalog Number: C11103

Concentration: 1mg/ml

Swiss-Prot No.: Q15743

Other Names: Sphingosylphosphorylcholine receptor, Ovarian cancer G-protein coupled receptor 1, OGR-1, G-protein coupled receptor 68, GPR12A, OGR1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human SPR1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: SPR1 antibody detects endogenous levels of total SPR1 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

An S., FEBS Lett. 375:121-124(1995).

Xu Y., Genomics 35:397-402(1996).

King M.M., Submitted (APR-2004) to the EMBL/GenBank/DDBJ databases.

