



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

CEBPE Antibody [#C11104]

Catalog Number: C11104

Concentration: 1mg/ml

Swiss-Prot No.: Q15744

Other Names: CCAAT/enhancer-binding protein epsilon, C/EBP epsilon, CEBPE

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human CEBPE.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CEBPE antibody detects endogenous levels of total CEBPE protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Antonson P., Genomics 35:30-38(1996).

Chumakov A.M., Mol. Cell. Biol. 17:1375-1386(1997).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

