



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

Smad1/5/9 Antibody [#C11110]

Catalog Number: C11110

Concentration: 1mg/ml

Swiss-Prot No.: Q15797

Other Names: BSP-1, BSP1, Dwarf-A, Dwf-A, JV4-1, MADH1, MADR1, Mad-related protein 1, Mad1, Mothers against DPP homolog 1, Mothers against decapentaplegic homolog 1, Mothers-against-DPP-related-1, SMAD 1, Smad 1, Transforming growth factor- beta signaling protein-1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against

synthesized peptide derived from N-terminal of human Smad1/5/9.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: Smad1/5/9 antibody detects endogenous levels of total Smad1/5/9 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Riggins G.J., Nat. Genet. 13:347-349(1996).

Liu F., Nature 381:620-623(1996).

Hoodless P.A., Cell 85:489-500(1996).

