



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ZNF187 Antibody [#C11138]

Catalog Number: C11138

Concentration: 1mg/ml

Swiss-Prot No.: Q16670

Other Names: SRE-ZBP; zinc finger and SCAN domain-containing 26; zinc finger protein 187; ZN187; ZSCAN26

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ZNF187.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ZNF187 antibody detects endogenous levels of total ZNF187 protein.

Reactivity: Human


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:5000

References:


Mungall A.J., Nature 425:805-811(2003).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Attar R.M., Mol. Cell. Biol. 12:2432-2443(1992).



Western blot analysis of extracts from Jurkat cells and HUVEC cells, using ZNF187 antibody.



Immunofluorescence analysis of A549 cells, using ZNF187 antibody.