



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

DGKD Antibody [#C11140]

Catalog Number: C11140

Concentration: 1mg/ml

Swiss-Prot No.: Q16760

Other Names: Diacylglycerol kinase delta, EC 2.7.1.107, Diglyceride kinase delta, DGK-delta, DAG kinase delta, 130 kDa diacylglycerol kinase, DGKD, KIAA0145

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

DGKD.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: DGKD antibody detects endogenous levels of total DGKD protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Sakane F., J. Biol. Chem. 277:43519-43526(2002).

Nagase T., DNA Res. 2:167-174(1995).

Sakane F., J. Biol. Chem. 271:8394-8401(1996).

