



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

AAK1 Antibody [#C11150]

Catalog Number: C11150

Concentration: 1mg/ml

Swiss-Prot No.: Q2M2I8

Other Names: AP2-associated protein kinase 1, EC 2.7.11.1, Adaptor-associated kinase 1, AAK1, KIAA1048

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human AAK1.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: AAK1 antibody detects endogenous levels of total AAK1 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Kikuno R., DNA Res. 6:197-205(1999).

Hillier L.W., Nature 434:724-731(2005).

Beausoleil S.A., Proc. Natl. Acad. Sci. U.S.A. 101:12130-12135(2004).

