



Web: [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
Order: [order@dijibio.com](mailto:order@dijibio.com)  
Tel: (86) 519 8805 0026

## ADCK5 Antibody [#C11156]

**Catalog Number:** C11156

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q3MIX3

**Other Names:** Uncharacterized aarF domain-containing protein kinase 5, EC 2.7.11.-, ADCK5

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human ADCK5.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** ADCK5 antibody detects endogenous levels of total ADCK5 protein.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Nusbaum C., Nature 439:331-335(2006).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

