



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

QORX Antibody [#C11167]

Catalog Number: C11167

Concentration: 1mg/ml

Swiss-Prot No.: Q53FA7

Other Names: Shugoshin-like 1, hSgo1, Serologically defined breast cancer antigen NY-BR-85, SGO1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human QORX.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: QORX antibody detects endogenous levels of total QORX protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

McGuinness B.E., PLoS Biol. 3:433-449(2005).

Kitajima T.S., Curr. Biol. 15:353-359(2005).

Tang Z., Proc. Natl. Acad. Sci. U.S.A. 101:18012-18017(2004).

