



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## SGOL1 Antibody [#C11176]

**Catalog Number:** C11176

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q5FBB7

**Other Names:** Shugoshin-like 1, hSgo1, Serologically defined breast cancer antigen NY-BR-85, SGO1

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human SGOL1.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** SGOL1 antibody detects endogenous levels of total SGOL1 protein.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

McGuinness B.E., PLoS Biol. 3:433-449(2005).

Kitajima T.S., Curr. Biol. 15:353-359(2005).

Tang Z., Proc. Natl. Acad. Sci. U.S.A. 101:18012-18017(2004).

