



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

MARK4 Antibody [#C11511]

Catalog Number: C11511

Concentration: 1mg/ml

Swiss-Prot No.: Q96L34

Other Names: MAP/microtubule affinity-regulating kinase 4, EC 2.7.11.1, MAP/microtubule affinity-regulating kinase-like 1, MARK4, KIAA1860, MARKL1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human MARK4.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MARK4 antibody detects endogenous levels of total MARK4 protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 100%, Positives = 100%)


Applications: IF: 1:100~1:500 ELISA: 1:40000

References:

Kato T., Neoplasia 3:4-9(2001).

Nagase T., DNA Res. 8:85-95(2001).

Ota T., Nat. Genet. 36:40-45(2004).



Peptide

-

+

Immunofluorescence analysis of A549 cells, using
MARK4 antibody.