



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TREF1 Antibody [#C11526]

Catalog Number: C11526

Concentration: 1mg/ml

Swiss-Prot No.: Q96PN7

Other Names: dJ139D8.5; DJ139D8.5.1; TREF1; Zinc finger transcription factor TReP-132; Zinc finger transcription regulating protein TReP-132, variant 1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human TREF1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TREF1 antibody detects endogenous levels of total TREF1 protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 100%, Positives = 100%)


Applications: IHC: 1:50~1:100 ELISA: 1:1000

References:

Gizard F., J. Biol. Chem. 276:33881-33892(2001).

Mungall A.J., Nature 425:805-811(2003).

Beausoleil S.A., Proc. Natl. Acad. Sci. U.S.A. 101:12130-12135(2004).



Peptide

-

+

Immunohistochemistry analysis of paraffin-embedded human brain tissue, using TREF1 antibody.