



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

LMTK3 Antibody [#C11528]

Catalog Number: C11528

Concentration: 1mg/ml

Swiss-Prot No.: Q96Q04

Other Names: Serine/threonine-protein kinase LMTK3
[Precursor], EC 2.7.11.1, Lemur tyrosine kinase 3, LMTK3,
KIAA1883, TYKLM3

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered
saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl,
0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against
synthesized peptide derived from internal of human

LMTK3.

Purification: The antibody was affinity-purified from
rabbit antiserum by affinity-chromatography using
epitope-specific immunogen.

Specificity: LMTK3 antibody detects endogenous levels
of total LMTK3 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:1000 ELISA: 1:20000

References:

Nagase T., DNA Res. 8:179-187(2001).

The MGC Project Team, Genome Res.
14:2121-2127(2004).

