



Web: [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
Order: [order@dijibio.com](mailto:order@dijibio.com)  
Tel: (86) 519 8805 0026

## AL2S7 Antibody [#C11529]

**Catalog Number:** C11529

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q96Q40

**Other Names:** A2S7, AL2S7, ALS2CR7, Amyotrophic lateral sclerosis 2 chromosomal region candidate gene protein 7, EC 2.7.11.22, Serine/threonine-protein kinase ALS2CR7

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human

AL2S7.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** AL2S7 antibody detects endogenous levels of total AL2S7 protein.

**Reactivity:** Human, Mouse

**Applications:** WB: 1:500~1:1000 ELISA: 1:5000

**References:**

Hadano S., Nat. Genet. 29:166-173(2001).

Ota T., Nat. Genet. 36:40-45(2004).

Hillier L.W., Nature 434:724-731(2005).

