



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

MAP3K3 Antibody [#C11574]

Catalog Number: C11574

Concentration: 1mg/ml

Swiss-Prot No.: Q99759

Other Names: EC 2.7.11.25, M3K3, MAP3K3, MAPK/ERK kinase kinase 3, MAPKKK3, MEK kinase 3, MEKK 3, Mitogen-activated protein kinase kinase kinase 3, kinase MEKK3

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

MAP3K3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MAP3K3 antibody detects endogenous levels of total MAP3K3 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Ellinger-Ziegelbauer H.C., J. Biol. Chem. 272:2668-2674(1997).
Nakamura K., J. Biol. Chem. 278:36989-36992(2003).
Huang Q., Nat. Immunol. 5:98-103(2004).

