



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

FOXB1/2 Antibody [#C11579]

Catalog Number: C11579

Concentration: 1mg/ml

Swiss-Prot No.: Q99853/Q5VYV0

Other Names: Forkhead box protein B1/2, Transcription factor FKH-5, FOXB1/2, FKH5

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human FOXB1/2.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: FOXB1/2 antibody detects endogenous levels of total FOXB1/2 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Bingle C.D., Submitted (JUN-1998) to the EMBL/GenBank/DDBJ databases.

Cascio S.M., Submitted (MAR-1998) to the EMBL/GenBank/DDBJ databases.

Bingle C.D., Submitted (JAN-1997) to the EMBL/GenBank/DDBJ databases.

