



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TINF2 Antibody [#C11607]

Catalog Number: C11607

Concentration: 1mg/ml

Swiss-Prot No.: Q9BSI4

Other Names: TRF1-interacting nuclear protein 2,
TERF1-interacting nuclear factor 2, TINF2, TIN2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human TINF2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TINF2 antibody detects endogenous levels of total TINF2 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Kim S.-H., Nat. Genet. 23:405-412(1999).

Ota T., Nat. Genet. 36:40-45(2004).

Loayza D., Nature 423:1013-1018(2003).

