



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

MMTA2 Antibody [#C11613]

Catalog Number: C11613

Concentration: 1mg/ml

Swiss-Prot No.: Q9BU76

Other Names: hMMTAG2, Multiple myeloma tumor-associated protein 2, MMTAG2, C1orf35

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human MMTA2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MMTA2 antibody detects endogenous levels of total MMTA2 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Tian J.Y., Sheng Wu Hua Xue Yu Sheng Wu Wu Li Xue Bao 35:143-148(2003).

Ota T., Nat. Genet. 36:40-45(2004).

Gregory S.G., Nature 441:315-321(2006).

