



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

MKKNK1 Antibody [#C11615]

Catalog Number: C11615

Concentration: 1mg/ml

Swiss-Prot No.: Q9BUB5

Other Names: MAP kinase-interacting serine/threonine-protein kinase 1, EC 2.7.11.1, MAP kinase signal-integrating kinase 1, Mnk1, MKKNK1, MNK1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

MKKNK1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MKKNK1 antibody detects endogenous levels of total MKKNK1 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Fukunaga R., EMBO J. 16:1921-1933(1997).

Knauf U., Mol. Cell. Biol. 21:5500-5511(2001).

O'Loghlen A., Exp. Cell Res. 299:343-355(2004).

