



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## UTP14A Antibody [#C11624]

**Catalog Number:** C11624

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q9BVJ6

**Other Names:** U3 small nucleolar ribonucleoprotein, homolog A; UT14A; UTP14

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human UTP14A.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** UTP14A antibody detects endogenous levels of total UTP14A protein.

**Reactivity:** Human, Mouse

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Ross M.T., Nature 434:325-337(2005).

Scherl A., Mol. Biol. Cell 13:4100-4109(2002).

Rohozinski J., Proc. Natl. Acad. Sci. U.S.A. 101:11695-11700(2004).

