



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

NADAP Antibody [#C11628]

Catalog Number: C11628

Concentration: 1mg/ml

Swiss-Prot No.: Q9BWU0

Other Names: Kidney anion exchanger adapter protein, Solute carrier family 4 anion exchanger member 1 adapter protein, Lung cancer oncogene 3 protein, HLC-3, SLC4A1AP, HLC3, Kanadaptin

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

NADAP.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: NADAP antibody detects endogenous levels of total NADAP protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Kittanakom S., Mol. Membr. Biol. 21:395-402(2004).

Ota T., Nat. Genet. 36:40-45(2004).

Hillier L.W., Nature 434:724-731(2005).

