



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

NARG1 Antibody [#C11637]

Catalog Number: C11637

Concentration: 1mg/ml

Swiss-Prot No.: Q9BXJ9

Other Names: N-terminal acetyltransferase, Protein tubedown-1, Tbdn100, Gastric cancer antigen Ga19

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human NARG1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: NARG1 antibody detects endogenous levels of total NARG1 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Line A., Br. J. Cancer 86:1824-1830(2002).

Willis D.M., J. Biol. Chem. 277:37280-37291(2002).

Fluge O., Oncogene 21:5056-5068(2002).

