



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## AN30A Antibody [#C11638]

**Catalog Number:** C11638

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q9BXX3

**Other Names:** Serologically defined breast cancer antigen NY-BR-1, Ankyrin repeat domain-containing protein 30A, ANKRD30A

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human

AN30A.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** AN30A antibody detects endogenous levels of total AN30A protein.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:1000

**References:**

Jaeger D., Cancer Res. 61:2055-2061(2001).

Deloukas P., Nature 429:375-381(2004).

Hahn Y., Gene 366:238-245(2006).

