



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

KPBB Antibody [#C11458]

Catalog Number: C11458

Concentration: 1mg/ml

Swiss-Prot No.: Q93100

Other Names: Phosphorylase b kinase regulatory subunit beta, Phosphorylase kinase subunit beta, PHKB

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human KPBB.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: KPBB antibody detects endogenous levels of total KPBB protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Wuellrich-Schmoll A., Eur. J. Biochem. 238:374-380(1996).

van den Berg I.E.T., Am. J. Hum. Genet. 61:539-546(1997).

