



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## CNKR1 Antibody [#C11460]

**Catalog Number:** C11460

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q969H4

**Other Names:** Connector enhancer of KSR1; hCNK1;  
Connector enhancer of KSR-like; CNK homolog protein 1

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human CNKR1.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** CNKR1 antibody detects endogenous levels of total CNKR1 protein.

**Reactivity:** Human, Mouse, Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:10000

**References:**

Therrien M., Cell 95:343-353(1998).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Jaffe A.B., Mol. Cell. Biol. 24:1736-1746(2004).

