



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

MKL1 Antibody [#C11466]

Catalog Number: C11466

Concentration: 1mg/ml

Swiss-Prot No.: Q969V6

Other Names: Myocardin-related transcription factor A,
MRTF-A, Megakaryoblastic leukemia 1 protein,
Megakaryocytic acute leukemia protein

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered
saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl,
0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against
synthesized peptide derived from N-terminal of human
MKL1.

Purification: The antibody was affinity-purified from
rabbit antiserum by affinity-chromatography using
epitope-specific immunogen.

Specificity: MKL1 antibody detects endogenous levels of
total MKL1 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Ma Z., Nat. Genet. 28:220-221(2001).

Mercher T., Proc. Natl. Acad. Sci. U.S.A.
98:5776-5779(2001).

Nagase T., DNA Res. 7:65-73(2000).

