



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TR19L Antibody [#C11468]

Catalog Number: C11468

Concentration: 1mg/ml

Swiss-Prot No.: Q969Z4

Other Names: Tumor necrosis factor receptor superfamily member 19L [Precursor], Receptor expressed in lymphoid tissues, RELT, TNFRSF19L

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human TR19L.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TR19L antibody detects endogenous levels of total TR19L protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Sica G.L., Blood 97:2702-2707(2001).

Ota T., Nat. Genet. 36:40-45(2004).

Cusick J.K., Biochem. Biophys. Res. Commun. 340:535-543(2006).

