



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TP53INP1 Antibody [#C11471]

Catalog Number: C11471

Concentration: 1mg/ml

Swiss-Prot No.: Q96A56

Other Names: Tumor protein p53-inducible nuclear protein 1, p53-dependent damage-inducible nuclear protein 1, p53DINP1, Stress-induced protein, TP53INP1, P53DINP1, SIP

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human

TP53INP1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TP53INP1 antibody detects endogenous levels of total TP53INP1 protein.

Reactivity: Human


Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Okamura S., Mol. Cell 8:85-94(2001).

Tomasini R., Eur. J. Cell Biol. 81:294-301(2002).

Tomasini R., J. Biol. Chem. 278:37722-37729(2003).



For Research Use Only