



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

IP3KC Antibody [#C11487]

Catalog Number: C11487

Concentration: 1mg/ml

Swiss-Prot No.: Q96DU7

Other Names: Inositol-trisphosphate 3-kinase C, EC 2.7.1.127, Inositol 1,4,5-trisphosphate 3-kinase C, InsP 3-kinase C, IP3K-C, ITPKC, IP3KC

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human IP3KC.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: IP3KC antibody detects endogenous levels of total IP3KC protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Dewaste V., Biochem. J. 352:343-351(2000).

Dewaste V., Biochem. J. 374:41-49(2003).

Structural genomics consortium (SGC), Submitted (OCT-2005) to the PDB data bank.

