



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

HES6 Antibody [#C11496]

Catalog Number: C11496

Concentration: 1mg/ml

Swiss-Prot No.: Q96HZ4

Other Names: Transcription cofactor HES-6, Hairy and enhancer of split 6, C-HAIRY1, HES6

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human HES6.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: HES6 antibody detects endogenous levels of total HES6.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Bae S.-K., Development 127:2933-2943(2000).

Vasiliauskas D., Mech. Dev. 98:133-137(2000).

Ota T., Nat. Genet. 36:40-45(2004).

