



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

BTBD6 Antibody [#C11508]

Catalog Number: C11508

Concentration: 1mg/ml

Swiss-Prot No.: Q96KE9

Other Names: Lens BTB domain protein, BTB/POZ domain-containing protein 6, BTBD6, BDPL

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from interf human BTBD6.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: BTBD6 antibody detects endogenous levels of total BTBD6 protein.


Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 100%, Positives = 100%)

Applications: IHC: 1:50~1:100 ELISA: 1:10000

References:

Wistow G., Mol. Vis. 8:171-184(2002).

The MGC Project Team, Genome Res. 14:2121-2127(2004).



Peptide

-

+

Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using BTBD6 antibody.