



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

FAM84B Antibody [#C11509]

Catalog Number: C11509

Concentration: 1mg/ml

Swiss-Prot No.: Q96KN1

Other Names: BCMP101; breast cancer membrane protein 101; FA84B; family with sequence similarity 84 member B; NSE2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human FAM84B.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: FAM84B antibody detects endogenous levels of total FAM84B protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Hughes P.J., Submitted (OCT-2001) to the EMBL/GenBank/DDBJ databases.

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Huang X.-P., Cancer Genet. Cytogenet. 165:20-24(2006).

