



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

MTLR Antibody [#G397]

Catalog Number: G397

Concentration: 1mg/ml

Swiss-Prot No.: O43193

Other Names: G protein-coupled receptor 38; G-protein coupled receptor 38; GPR38; MLNR; Motilin receptor; MTLR; MTLR1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human MTLR.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MTLR antibody detects endogenous levels of total MTLR protein.

Reactivity: Human


Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

McKee K.K., Genomics 46:426-434(1997).

Feighner S.D., Science 284:2184-2188(1999).

Dunham A., Nature 428:522-528(2004).



For Research Use Only