



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

MTR1B Antibody [#G399]

Catalog Number: G399

Concentration: 1mg/ml

Swiss-Prot No.: P49286

Other Names: Melatonin receptor type 1B; Mel1b melatonin receptor; Mel-1B-R; MTNR1B

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human MTR1B.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MTR1B antibody detects endogenous levels of total MTR1B protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Reppert S.M., Proc. Natl. Acad. Sci. U.S.A. 92:8734-8738(1995).

Ebisawa T., Neurosci. Lett. 280:29-32(2000).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

