



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

NK1R Antibody [#G402]

Catalog Number: G402

Concentration: 1mg/ml

Swiss-Prot No.: P25103

Other Names: Substance-P receptor; SPR; NK-1 receptor; NK-1R; Tachykinin receptor 1; TACR1; NK1R, TAC1R

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human NK1R.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: NK1R antibody detects endogenous levels of total NK1R protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Takeda Y., Biochem. Biophys. Res. Commun. 179:1232-1240(1991).

Hopkins B., Biochem. Biophys. Res. Commun. 180:1110-1117(1991).

Gerard N.P., Biochemistry 30:10640-10646(1991).

