



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** +86 519 8805 0026

## NPY1R Antibody [#G409]

**Catalog Number:** G409

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P25929

**Other Names:** Neuropeptide Y receptor type 1; NPY1-R;  
NPY1R; NPYR, NPYY1

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human NPY1R.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** NPY1R antibody detects endogenous levels of total NPY1R protein.

**Reactivity:** Human, Mouse, Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:1000

**References:**

Larhammar D., J. Biol. Chem. 267:10935-10938(1992).

Herzog H., Proc. Natl. Acad. Sci. U.S.A. 89:5794-5798(1992).

Herzog H., J. Biol. Chem. 268:6703-6707(1993).

