



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

AKR1CL1 Antibody [#C14397]

Catalog Number: C14397

Concentration: 1mg/ml

Swiss-Prot No.: Q5T2L2

Other Names: Aldo-keto reductase family 1 member C-like protein 1; EC 1.1.1.-

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human AKR1CL1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: AKR1CL1 antibody detects endogenous levels of total AKR1CL1 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)

Applications: IF: 1:100~1:500 ELISA: 1:40000

References:

Ota T., Nat. Genet. 36:40-45(2004).

Deloukas P., Nature 429:375-381(2004).

The MGC Project Team;

Genome Res. 14:2121-2127(2004).

