



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

AKR1CL2 Antibody [#C14398]

Catalog Number: C14398

Concentration: 1mg/ml

Swiss-Prot No.: Q96JD6

Other Names: akcl2; akr1cl2; akr1cl2 protein; akrdc1; aldo-keto reductase family 1, member c-like 2; aldo-keto reductase loopadr; aldo-keto reductase related protein 1; aldo-keto reductase related protein 2; aldo-keto reductase related protein 3; aldo-keto reductase related protein 4; loopadr; mgc10612

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

AKR1CL2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: AKR1CL2 antibody detects endogenous levels of total AKR1CL2 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Nishinaka T., Chem. Biol. Interact. 143:299-305(2003).

Hyndman D.J., Submitted (MAY-2000) to the EMBL/GenBank/DDBJ databases.

The MGC Project Team; Genome Res. 14:2121-2127(2004).

