



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

AKR1B1 Antibody [#C14400]

Catalog Number: C14400

Concentration: 1mg/ml

Swiss-Prot No.: P15121

Other Names: ADR; aldehyde reductase 1; aldo-keto reductase family 1, member B1; aldose reductase; ALDR; ALDR1; AR; EC 1.1.1.21; Lii5-2 CTCL tumor antigen; low Km aldose reductase

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

AKR1B1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: AKR1B1 antibody detects endogenous levels of total AKR1B1 protein.

Reactivity: Human, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Bohren K.M., J. Biol. Chem. 264:9547-9551(1989).

Chung S., J. Biol. Chem. 264:14775-14777(1989).

Graham A., Nucleic Acids Res. 17:8368-8368(1989).

