



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

AMPD1 Antibody [#C14443]

Catalog Number: C14443

Concentration: 1mg/ml

Swiss-Prot No.: P23109

Other Names: AMP deaminase 1; EC 3.5.4.6;

Myoadenylate deaminase; AMP deaminase isoform M

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human AMPD1.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: AMPD1 antibody detects endogenous levels of total AMPD1 protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 100%, Positives = 100%);
Rat (Identities = 100%, Positives = 100%)

Applications: IHC: 1:50~1:100 ELISA: 1:40000

References:

Gregory S.G., Nature 441:315-321(2006).
Morisaki H., Hum. Mutat. 16:467-472(2000).
Sjoeblom T., Science 314:268-274(2006).

