



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

AMPH Antibody [#C14446]

Catalog Number: C14446

Concentration: 1mg/ml

Swiss-Prot No.: P49418

Other Names: AMPH; AMPH1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human AMPH.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: AMPH antibody detects endogenous levels of total AMPH protein.

Reactivity: Human (Identities = 100%, Positives = 100%);

Mouse (Identities = 100%, Positives = 100%);

Rat (Identities = 100%, Positives = 100%)

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

David C., FEBS Lett. 351:73-79(1994).

Yamamoto R., Hum. Mol. Genet. 4:265-268(1995).

Floyd S.R., Mol. Med. 4:29-39(1998).

