



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

APLP2 Antibody [#C14458]

Catalog Number: C14458

Concentration: 1mg/ml

Swiss-Prot No.: Q06481

Other Names: Amyloid protein homolog; Amyloid-like protein 2 precursor; APLP2; APP2; APPH; APPL2; CDEBP; CDEI-box binding protein

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human APLP2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: APLP2 antibody detects endogenous levels of total APLP2 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)
Mouse (Identities = 100%, Positives = 100%);
Rat (Identities = 100%, Positives = 100%)

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Sprecher C.A., Biochemistry 32:4481-4486(1993).

Wasco W., Nat. Genet. 5:95-99(1993).

von der Kammer H., DNA Cell Biol. 13:1137-1143(1994).

