



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

AIG1 Antibody [#C14468]

Catalog Number: C14468

Concentration: 1mg/ml

Swiss-Prot No.: Q9NVV5

Other Names: Androgen-induced protein 1; AIG-1;
CGI-103

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human AIG1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: AIG1 antibody detects endogenous levels of total AIG1 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Seo J., Mol. Cells 11:35-40(2001).

Lai C.-H., Genome Res. 10:703-713(2000).

Ota T., Nat. Genet. 36:40-45(2004).

