



Web: [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
Order: [order@dijibio.com](mailto:order@dijibio.com)  
Tel: (86) 519 8805 0026

## ANGPTL7 Antibody [#C14484]

**Catalog Number:** C14484

**Concentration:** 1mg/ml

**Swiss-Prot No.:** O43827

**Other Names:** angiopoietin-like 7; ANGL7; AngX; CDT6; cornea-derived 6

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from C-terminal of human ANGPTL7.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** ANGPTL7 antibody detects endogenous levels of total ANGPTL7 protein.

**Reactivity:** Human, Mouse

**Applications:** WB: 1:500~1:3000 ELISA: 1:10000

**References:**

Peek R., Invest. Ophthalmol. Vis. Sci. 39:1782-1788(1998).

Stover C., Genes Immun. 2:119-127(2001).

Clark H.F., Genome Res. 13:2265-2270(2003).

